

Setup New Session inside Motive

Before you start capturing you need to setup your session.

first you have to open the data-pane:



then press the plus-icon to select a new session folder



you need to create a new session inside your projetct on your filer:

I:\Innolab\03_projects\180101_Pilot_DemoProject\01_capturing

if necessary create a new folder with a distinctive name with a date,

Session_180329



then you need to remove all the other sessions from previous uses:



ATTENTION: PLEASE make sure you only REMOVE the folders. DO NOT DELETE THEM!!!



Set the blue flag to make sure your new session folder is active and all new takes will be stored inside this folder:



For more detailed information you can check the [official documentation](#):

From:

<https://wiki.zhdk.ch/IASpace/> - **immersive art space**

Permanent link:

https://wiki.zhdk.ch/IASpace/doku.php?id=motive_sessions&rev=1522155127

Last update: **2018/03/27 14:52**

