2025/07/22 02:11 1/1 Setup Tracking Suit

## Setup the tracking suit

Use the tracking suits ONLY with the black cotton-bodies underneath. It is the users responsibility to wash the bodies and bring them back in a timely manner.

Choose a suitable size suit for the actor. Make sure the suit fits tight and the top and pants cannot shift. Take off all the markers that are still on the suit.

We are using the Markerset "BaseLine + Hinged Toe" with 41 markers:





feet: 5 (x 2)legs: 3 (x 2)

• hips: 4 (2 front - 2 back)

back: 3 (triangle)shoulder: 2 (x 2)

chest: 1arm: 2 (x 2)hands: 3 (x 2)head: 3

take care at the joints. set the markers in such a way, that the marker defines the axis of the joint and doesn't move when the joint is used. → ankle, knee, shoulder, elbow.

let the actor stand in the center of the space and take the T-pose.

the system should detect now all the 41 markers.

give it a name and press create.

If everthing went right, the virtual character should appear inside perpective view. make sure the settings allow the display of the skeletons:



From:

https://wiki.zhdk.ch/IASpace/ - immersive art space

Permanent link:

https://wiki.zhdk.ch/IASpace/doku.php?id=motive\_trackingsuit&rev=1522161148

Last update: 2018/03/27 16:32

